

**Silicon Beads Lab Test Report
26 June 2007**

Sample Identity- Canty Sample ID 1340

- 1 Large Product
- 2 Small Product

Purpose:

This report details the particle size testing done in the Canty lab in order to ascertain size and repeatability of measurement.

Lab Setup:



Canty SolidSizer system was used for these dry products.
Datasheet Link: [SolidSizer](#)
Video Link: [Video](#)

Typical Image:



Results:

**Sample 1340-1: Large Particle Product
Calibration – 77.4 microns per pixel**



Sample 2: Small Particle Product
Calibration: 31.35 microns per pixel

